

Molift Partner 255

The flagship of Molift's mobile hoists - with vertical hoisting pattern for excellent stability

Ideal solution for critical care

Molift Partner 255 provides an ideal solution for hoisting in acute and critical care areas in both hospitals and care facilities. The straight vertical column in combination with the fixed sling bar provides excellent stability. It maintains the user at the same distance from the column throughout the hoisting. With an outstanding hoisting range, Molift Partner 255 allows hoisting from the floor as well as high surfaces.

Molift Partner 255 is one of the lightest in its category and combined with excellent construction, ergonomic design and the caster solution, it becomes very easy to manoeuvre even with a lifting capacity of 255 kg.

Molift Partner 255 is available in several variants for the possibility of optimizing the product according to needs; low base chassis, supporting arms for gait/walk training and weight scale.



Features

- High safe working load - 255 kg
- Fixed sling bar and straight vertical column for safe and stable hoisting
- Low weight and easy to manoeuvre
- Suitable for use with Molift stretcher

| Product | Item no. | Dimensions | Technical |
|--|----------|--|--|
| Molift Partner 255 incl. sling bar | 13000 | Length 1285 mm | Battery NiMh 26,4 V-2,6 Ah |
| For complete variant matrix, contact your local distributor. | | Chassis width (outer) 700 mm | Battery charge time approx. 3 hours |
| Accessories | Item no. | Chassis width (inner) 550 mm | Hoisting speed 30 mm/sec |
| Casters ø100mm with directional lock | 1320140 | Chassis width, open legs (inner) 1065 mm | Rated performance 50 hoists/75 kg/500 mm |
| Weight scale (see page for sling bars) | | Hoisting height 390-1680 mm | Material |
| Sling bars (see page for sling bars) | | Hoisting interval 1290 mm | Aluminium, Steel |
| Stretchers (see page for stretchers) | | Leg height 115 mm | |
| For more information visit etac.com | | Weight | |
| | | Total weight 42.5 kg | |
| | | | |